

CERVICAL INTERBODY CAGE PM CAGE® INSTRUCTIONS FOR USE

Manufacturer: NEURO FRANCE Implants
ZA LE BOURG – 25 RUE DES ECOLES
41160 LA VILLE AUX CLERCS – FRANCE
Telephone: (33) 02 54 80 90 90
Fax: (33) 02 54 80 83 33



Important information for the medical practitioner:

The practicing physician should carefully read the instructions and advice as well as the recommendations before using the cervical interbody cage "PM Cage®".

1. Description of the system

- The system is composed of PM cage® (Small or large).
- The PM Cage® system is made of either titanium alloy Ti 6-Al 4-V according to the ISO 5832-3 standard, or of PEEK Optima® (with radiological markers made of titanium alloy Ti 6-Al 4-V).

2. Indications

The cervical interbody cage "PM Cage®" system is for the treatment of spine diseases such as:

- Cervicobrachial neuralgia (C.B.N.)
- Degenerative spine diseases
- Recurrent spinal disc herniation
- Cervical instability
- Spinal trauma
- Uncarthrosis.

3. Contraindications

The choice of the treatment is determined by the practitioner who has the required skills, knowledge and experience.

Factors which may compromise the implantation include (but are not limited to):

- Malignant vertebral diseases
- Infection,
- Allergy to components of the titanium alloy Ti 6-Al 4-V and to PEEK Optima®
- Impatient patient overweight
- Hyperactivity/mental disorders
- Severe vertebral osteoporosis, osteomyelitis
- Pregnancy.

4. Use

- Please consult the surgical technique before use.
- The surgical approach is anterior.

5. Warning / Operating precautions

- The surgeon must be fully experienced in the use of the PM Cage® implants, the application method, the instruments and the surgical technique.
- The use of similar implants in conjunction with the PM Cage® system is proscribed.
- The use of other instruments and/or the incorrect use of the instruments required to fit the PM Cage® implants may lead to a deterioration of the implants that will limit the optimization of the desired function. It is therefore proscribed.

- Although allergic reactions are extremely rare, we recommend to make an allergy check-up on patients with a particularly strong predisposition to allergies.
- The PM Cage® system must be used under perfectly sterile conditions.
- The implants are for single use only and should never be reused. Reuse may lead to the degradation of device performances, contaminations and cross-infections.
- All implants extracted from the patient, damaged or having been incorrectly used should be disposed of once it has been soiled or contaminated with blood, body tissues.
- There are no known risks of mutual interferences between the implant and other medical devices.
- The patient should ask for medical advice before entering potentially dangerous environments (such as electromagnetic or magnetic fields including magnetic resonance environments), which may affect the performances of the implant.
- The implants are made of titanium alloy Ti 6-Al 4-V or PEEK Optima® which are non-magnetic materials. Patients can be exposed to electromagnetic or magnetic fields without any real risks. The patient must notify that he has implant before any exposure to electromagnetic or magnetic fields.
- The patient has to be informed of the risks of the surgical intervention by the medical staff. The patient should follow the advice and recommendations from the surgeon (radiological examinations) during the postoperative phase.
- It is recommended to limit physical activity after the surgery to allow bone fusion. If not, the implant may break or otherwise be damaged necessitating revision surgery. The implant should not be exposed to extreme movements of mechanical vibrations for example.
- All modifications (aspect, pain...) at the implant site must be reported to the medical practitioner. The medical practitioner should also be informed of all type of incidents like a fall for instance, even when there are no visible signs at the implant site.

6. Adverse effects and complications

The potential adverse effects and complications are:

- Severe pseudarthrosis
- Surgical site infection
- Displacement or expulsion of the implant which may require revision surgery
- Neurological damage, breach of the dura mater, nerve root injury
- Vascular damage caused by the surgical procedure
- Stress fracture of the implant
- Death.

7. Storage Conditions

Implants must be stored in their original packaging in a dry environment and at room temperature.

8. Sterilization

The PM cage® implants and the instrument set are supplied non-sterile. They must be sterilized under the responsibility of the person responsible for sterilization in the healthcare facility. Sterilization must be carried out in an autoclave.

Method: Sterilization by steam (moist heat)

Minimum duration: 18 minutes

Minimum temperature: 134°C

It is recommended to check the sterilization indicators and the integrity of the packaging before opening.

9. Disposal of implants

Implants must be disposed of according to the relevant hygiene and waste disposal guidelines of the healthcare facility for medical waste.

10. Additional Information

Put on the market: March 2005

Additional information can be obtained from NEURO FRANCE Implants.

Examples of symbols used on the labels:

	Batch number		Non sterile
	Catalogue number		Single use
	Manufacturer		CE mark with Notified Body identification number

Claims:

All claims and/or malfunction reports should be addressed immediately by telephone, fax or letter to NEURO FRANCE Implants.

Revised 01/02/2018