ThoRib® & TRIONYX®

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ThoRib® & TRIONYX®

INSTRUCTIONS FOR USE

Manufacturer: NEURO FRANCE Implants (NFI)

ZA LE BOURG – 25 RUE DES ECOLES
41160 LA VILLE AUX CLERCS – FRANCE
Telephone: (33) 02 54 80 90 90
Fax: (33) 02 54 80 83 33
Email: sav@neuro-france.net

Important information for the practitioner:

The practitioner must carefully read the instructions and guidelines, as well as the recommendations (presented in the surgical technique and in these instructions) before using the thoracic osteosynthesis system.

1. **Description of the materials**

- These instructions for use cover the generic group of thoracic osteosynthesis systems. They include the ThoRib® system, the ThoRib® PECTUS system and the TRIONYX® system.
- The systems are made up of:

Medical device	Intended purpose	Material
ThoRib®		
Costal staples	Solid fixation on a rib	Titanium T40 (standard ISO 5832-2)
Simple/double costal staples	Fracture reduction	
Sternal plates with notches	Sternal bridging on ribs	
Rods	Bridging one or several ribs	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
Connectors	Connection of a rib to a spinal rod	
Rod-staple junction screw with brake insert	Fixation of a rod or plate to a staple	Ti-6Al-4V titanium alloy (standard ISO 5832-3) PEEK Optima™ for the brake insert
ThoRib® PECTUS		
Pectus staples	Solid fixation on a rib	Titanium T40 (standard ISO 5832-2)
Pectus rods	Correction of sternal deformity	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
Rod-staple junction screw with long brake insert	Fixation of a rod to a staple	Ti-6Al-4V titanium alloy (standard ISO 5832-3) PEEK Optima™ for the brake insert
TRIONYX®		
Sternal staples	Solid fixation on a rib	Titanium T40 (standard ISO 5832-2)
Sternal plates	Partial or total replacement of the sternum	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
Rod-staple junction screw with long brake insert	Fixation of a plate to a staple	Ti-6Al-4V titanium alloy (standard ISO 5832-3) PEEK Optima™ for the brake insert

- These systems are available in non-sterile condition.
- The batch number is laser-printed on the implant. It is also printed on the product packaging.
- These implants are not intended to transport or introduce any substance into the human body.
- These implants do not contain any medicinal product, derivative of human blood or plasma, or tissue or cells of human or animal origin.
- Appropriate instruments are provided for implanting these systems. These ancillary products are reusable surgical instruments.
- These instruments are made of stainless steel or titanium. The handles are made of stainless steel, except for some which are made of silicone.
- They are delivered non-sterile. *Specific instructions for use are available for these instruments.*

2. **Intended purpose**

Medical device	Intended purpose
ThoRib®	Thoracic reconstruction after rib removal <i>In the case of trauma:</i> Creation of a costal bridge and synthesis of rib fractures.
ThoRib® PECTUS	Reduction of deformities and malformations of the chest wall (pectus excavatum, carinatum and arcuatum)
TRIONYX®	Massive thoracic reconstruction after removal of the sternum and/or the manubrium.

Under no circumstances should ThoRib® implants be used to treat deformities and malformations of the chest wall. Only the PECTUS range is designed for this type of indication.

3. **Expected clinical performance**

Medical device	Expected clinical performance
ThoRib®	Restoration of thoracic rigidity in over 90% of cases, following: - Wide resection of a chest wall bone tumour - Thoracic trauma
ThoRib® PECTUS	Correction of sternal deformity in over 90% of cases.
TRIONYX®	Restoration of thoracic rigidity in over 90% of cases, limiting post-operative morbidity and mortality.

4. **Expected clinical benefits**

Medical device	Expected clinical benefits
ThoRib®	Reduction in current clinical pain scores (VAS) of at least 80% at 12 months postoperatively.
ThoRib® PECTUS	Reduction in current clinical pain scores (VAS) of at least 80% at 12 months postoperatively.
TRIONYX®	Preservation of respiratory mechanics in over 90% of cases. Prevention of paradoxical respiration in over 90% of cases.

5. Indications

Medical device	Indications	Patient profile
ThoRib®	Severe thoracic trauma Chest wall and/or sternal tumour	Any adult patient (male or female) with one of the indications previously described and with no contraindications.
ThoRib® PECTUS	Deformities and malformations of the chest wall (pectus excavatum, carinatum and arcuatum)	
TRIONYX®	Sternal and/or manubrial bone tumour	

6. Contraindications

The following conditions may compromise the implantation (the list is not exhaustive):

- Severe osteoporosis
- Osteomalacia
- Brittle bone disease
- Active infectious process or major risk of infection
- Weakened immune system
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Hyperactivity
- Mental illness
- Allergy or intolerance to the implant materials
- Any case not described in the indications

7. Use

- Consult the surgical technique before use.

Medical device	Surgical approach
ThoRib®	Anterior and/or anterolateral approach
ThoRib® PECTUS	
TRIONYX®	

- The choice of staple size and/or TRIONYX® plate must take into account the size of the patient's bone resection in order to use implants of suitable lengths.

8. User profile

Due to the technical difficulties involved and the potential risk to the patient, the procedure should only be carried out by a surgeon who has been duly trained in thoracic surgery.

Training for surgeons in the placement of the ThoRib®, ThoRib® PECTUS and TRIONYX® systems is available. It shows how to handle surgical instruments and model implants. Reading the documentation and surgical techniques allows the surgeon to assimilate and understand the surgical risks associated with placement of the implant.

Sterilisation staff must refer to these instructions for use in particular to understand the cleaning and sterilisation procedures to be applied.

9. Warning/Surgical precautions

- The surgeon must be completely familiar with the implants, instruments and the surgical technique.
- The surgeon must take into account all of the factors likely to influence the performance of the device: number of ribs to be treated, weight, activity level, concomitant disorders, etc.
- Although allergic reactions to titanium are very rare, it is recommended that patients with particularly strong allergies be checked beforehand for any allergies of this type.
- The ThoRib®, ThoRib® PECTUS and TRIONYX® systems must be used under perfectly sterile surgical conditions.
- The components of the ThoRib®, ThoRib® PECTUS and TRIONYX® systems must never be combined with components of other devices.
- For the TRIONYX implant, a minimum of 4 staples (2 on each side) must be placed on the ribs for reasons of stability and sufficient mechanical resistance.
- Incorrect use of the instruments during implantation may lead to deterioration of the implants and therefore impact the expected performance.
- The instruments provided are designed to fit perfectly with the ThoRib®, ThoRib® PECTUS and TRIONYX® systems. It is therefore strongly contraindicated to use any other instruments. This may incur risks for patients (see Sec. 10) and damage the implants.
- Our implants are for single use only and must never be reused. Reuse may lead to reduced performance of the device, contamination and cross-infection.
- The implant must be disposed of if it has come into contact with blood, biological fluids or bodily tissues.
- The implant must be disposed of if it has been removed from the patient, damaged or if it has been used incorrectly.
- The implants are manufactured using Ti-6Al-4V titanium alloy, T40 titanium and PEEK Optima®, which are non-magnetic materials and have a geometry that is not likely to generate significant induced currents. Moreover, as they are fixed to the bone and/or tissue, the implants are unlikely to be moved. Their safety, specifically in terms of heating and implant migration, has been assessed on the basis of bibliographical data comparing the available data on devices of similar composition, shape and use. This assessment concluded that they were safe for use in patients undergoing MRI examinations between 1.5 and 3 Tesla.
- As a precaution, it is advisable to inform the MRI operator of the presence of these implants before these examinations.
- Titanium implants can generate metal-related artefacts, which must be taken into account when performing and interpreting imaging examinations.
- The surgeon must be notified of any change (appearance, pain, etc.) at the site of the implant.
- Any type of accident must be reported to the surgeon, even if there are no visible exterior indications at the site of the implant.

10. Adverse effects/complications

The potential adverse effects or complications are:

- Superficial or deep implantation site infection, skin necrosis
- Hypersensitivity reaction to the implants
- Pain, discomfort, stiffness
- Vascular, cardiac or pulmonary damage linked to the surgical procedure
- Displacement, fracture and/or expulsion of the implant requiring a repeat intervention
- Fatigue fracture of the material and/or bone
- Death

11. Storage conditions for implants

The implants must be stored in a dry place at room temperature.

12. Cleaning

The implants must be cleaned and disinfected before being sterilised. The purpose of this step is to reduce the number of viable micro-organisms on the implants prior to sterilisation.

The chemicals used during cleaning and disinfection must be compatible with the implant. The implants must not come into contact with products containing chlorine, phosphorous, formalin, fluorine, bleaching agents, fatty-acid detergents or any strong acid or alkaline solutions. To guarantee effective cleaning and disinfection, please follow the cleaning product instructions.

The following method must be used for decontamination and cleaning:

Equipment required: Washer-disinfector, enzymatic cleaner or detergent with a pH between 8 and 10, such as Dentalzym5.

Programme validated with the SMEG WD4060D washer-disinfector.

1. Remove the implants from their packaging and place them in a basket on their respective racks.
2. Place the basket in the washer-disinfector using the validated process so as to avoid dead zones (Caution: do not stack multiple baskets on top of each another. Cleaning in a washer-disinfector must generally be carried out without a cover).
3. Follow the programme below:

SMEG WD4060D washer	Duration (in minutes)	Temperature	Cleaning agent	Quantity
Pre-wash	3	Unheated	Tap water	-
Wash cycle 1	5	35°C	Dentalzym5	6 mL/litre of water (Volume of water: 6 L)
Rinse cycle	-	Unheated	Tap water	-
Wash cycle 2	5	75°C	Dentalzym5	6 mL/litre of water (Volume of water: 6 L)
Rinse cycle	-	Unheated	Tap water	-
Neutralisation	3	-	Dental AC5	10 mL/litre of water (Volume of water: 9 L)
Rinse cycle	3	-	Tap water	-
Thermal disinfection	5	90°C	Reverse osmosis water	-
Drying cycle	25	110°C	-	-

4. When the programme ends, remove the basket from the washer-disinfector and visually confirm that the wash has been completed.
5. Place the basket in folded packaging (envelope or Pasteur).

13. Sterilisation

As the implants are delivered non-sterile, the health institution is responsible for sterilising them.

Specialised containers are available to staff for storing the implants.

Sterilisation must be carried out in an autoclave, using one of the following protocols.

Method 1	Method 2
Steam sterilisation (moist heat) Recommended time: 18 minutes (under saturated steam pressure) Recommended temperature: 134°C	Steam sterilisation (moist heat) Recommended time: 4 minutes (under saturated steam pressure) Recommended temperature: 132°C

It is recommended that operating theatre staff check the integrity of packaging after sterilisation in the autoclave.

14. Disposing of implants

The removal and disposal of implants follow the same procedures used in the health institution as for hospital waste.

15. Medical device vigilance

Any complaint or notification of malfunction must be made immediately either by telephone, fax, email or post to NEURO FRANCE Implants.

NEURO FRANCE Implants and the competent authority of the Member State where the user and/or the patient is based must be notified of any serious incident involving the device.

16. Lifetime

The lifetime of this system is 18 years. This consists of:

- A maximum shelf life of 10 years
- And an implantation period of 8 years.

Please note that this lifetime is set based on mechanical, biological and clinical validations carried out on these medical devices.

17. Patient information

The surgeon must warn the patient of the risks of surgery and inform them of the potential adverse effects.

The patient must be informed that the devices cannot reproduce the flexibility, strength, reliability or durability of normal healthy bone. Therefore, if the patient's profession or activities subject the implant to excessive stresses (e.g. lifting loads or making major muscular efforts, exposure to excessive mechanical vibration movements, etc.), the implant may break or be damaged, which may require another operation. The surgeon must therefore discuss all of the physical and psychological limitations inherent in the use of these devices with the patient. The surgeon must also explain the need for regular medical follow-up during the post-operative phase.

Smoking has been shown to impede patients' recovery after surgery and affect bone healing and repair. Surgeons must inform patients of this fact and warn them of the potential consequences.

Patients should be aware that they must declare that they have an implant before any exposure to electromagnetic and magnetic environments.

For further information, a patient information letter is available on the NEURO FRANCE Implants website.

Health institutions and/or practitioners are asked to complete the patient card provided by NEURO FRANCE Implants. This card is intended to inform the patient of the type and lot number of the implanted medical device, the health institution and the surgeon.

18. Removal of implants

PECTUS implants are temporary internal fixation devices intended to stabilise the surgical site during the healing process.

After fusion, these implants no longer have any functional use and can be removed. The ancillary equipment supplied by NEURO FRANCE Implants can be used to remove implants.

Each decision to remove the device must take into account factors such as the risks of this surgical procedure for the patient, as well as the difficulty of removing the material, etc.

Removal may also be necessary in other cases, such as:

- Another operation for pseudarthrosis or other lack of bone fusion
- Migration of the implant, accompanied by pain and/or lesions
- Pain or abnormal sensations due to the presence of implants
- Infection or inflammatory reactions, etc.

19. Additional information

Medical device	Date placed on the market
ThoRib®	June 2012
ThoRib® PECTUS	June 2012
TRIONYX®	June 2016

Information on the clinical results of these systems can be found in the "Summary of Safety and Clinical Performance". This report is available in the European Medical Device Database (Eudamed)*, where it is linked to the Basic UDI-DI of different thoracic osteosynthesis systems:

Medical device	Basic UDI-DI system
ThoRib®	37005271THOSYSTEMRZ
ThoRib® PECTUS	37005271PECSYSTEMEY
TRIONYX®	37005271TRISYSTEMUR

The public Eudamed website is accessible via the following link: <https://ec.europa.eu/tools/eudamed>

Further information can be obtained from NEURO FRANCE Implants.

Explanation of the symbols used on the product labelling:

	Catalogue number		Protective packaging
	Batch code		Manufacturer and date of manufacture
	Do not re-use		Consult instructions for use
	Non-sterile		CE marking with the identification number of the notified body
	Medical device		Unique Device Identifier

These implants are Class IIb medical devices and bear the CE marking. Our notified body is DQS.

**Once the notified body has validated the report and the Eudamed site is fully operational.*

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Medical device	Intended purpose	Material
ThoRib®		
Costal staples	Solid fixation on a rib	Titanium T40 (standard ISO 5832-2)
Simple/double costal staples	Fracture reduction	
Sternal plates with notches	Sternal bridging on ribs	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
Rods	Bridging one or several ribs	
Connectors	Connection of a rib to a spinal rod	
Rod-staple junction screw with brake insert	Fixation of a rod or plate to a staple	Ti-6Al-4V titanium alloy (standard ISO 5832-3) PEEK Optima™ for the brake insert
ThoRib® PECTUS		
Pectus staples	Solid fixation on a rib	Titanium T40 (standard ISO 5832-2)
Pectus rods	Correction of sternal deformity	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
Rod-staple junction screw with long brake insert	Fixation of a rod to a staple	Ti-6Al-4V titanium alloy (standard ISO 5832-3) PEEK Optima™ for the brake insert
TRIONYX®		
Sternal staples	Solid fixation on a rib	Titanium T40 (standard ISO 5832-2)
Sternal plates	Partial or total replacement of the sternum	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
Rod-staple junction screw with long brake insert	Fixation of a plate to a staple	Ti-6Al-4V titanium alloy (standard ISO 5832-3) PEEK Optima™ for the brake insert

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17. Patient information

The surgeon must warn patients of the risks of surgery and inform them of the potential adverse effects.

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Health institutions and/or practitioners are asked to complete the implant card provided by NEURO FRANCE Implants. This card is intended to inform the patient of the type and lot number of the implanted medical device, the health institution and the surgeon.

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