CAGES STERILE

En - LUMBAR & CERVICAL INTERBODY CAGES (STERILE)

INSTRUCTIONS FOR USEp.2

Fr - CAGES INTERSOMATIQUES LOMBAIRES & CERVICALES (STERILE)

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Important information for the practitioner:

The practitioner must carefully read the instructions and guidelines, as well as the recommendations (presented in the surgical technique and in these instructions) before using a lumbar or cervical interbody cage.

1. Description of the materials

- These instructions for use cover the interbody cages generic group. It includes lumbar interbody cages supplied sterile (SHARK and PDP) and cervical interbody cages supplied sterile (PM CAGE®).
- These devices are made from PEEK Optima™ (with radiological markers made of Ti-6Al-4V titanium alloy).
- The batch number is marked on the implant. It is also printed on the product packaging.
- These implants are not intended to transport or introduce any substance into the human body.
- These implants do not contain any medicinal product, derivative of human blood or plasma, or tissue or cells of human or animal origin. Only the PE packaging may contain a small amount of animal fat derivatives, though the production of these derivatives complies with Regulation (EC) No. 999/2001. The risks that could result from the use of these additives are controlled and negligible.
- Appropriate instruments are provided for implanting these devices. These ancillary products are reusable surgical instruments.
- These instruments are made of stainless steel or titanium. The handles are made of stainless steel, except for some which are made of silicone.
- They are delivered non-sterile. *Specific instructions for use are available for these instruments.*

2. Intended purpose

These cages are designed for interbody fusion. These implants restore the height of the intervertebral space.

Note that lumbar interbody cages were not designed as "independent" implants, so they require the use of posterior instrumentation (such as pedicle screws or plates)

3. Expected clinical performance

- Fusion rate greater than 85% at 12 months postoperatively.
- Increase in disc height of at least 15%.

4. Expected clinical benefits

- Reduction in current clinical pain score (VAS) of at least 60% at 12 months postoperatively.
- Reduction in current clinical disability score (NDI/ODI) of at least 40% at 12 months postoperatively.

5. Indications

Medical device	Indications	Patient profile	Levels
SHARK	Degenerative disc disease Spondylolisthesis	Any adult patient (male or female) with one of the indications previously described and no contraindications.	L1 - S1
PDP			
PM CAGE®	Degenerative disc disease Cervico-brachial neuralgia (myelopathy - radiculopathy)		C2 – T1

6. Contraindications

The following conditions may compromise the implantation (the list is not exhaustive):

- Severe osteoporosis;
- Osteomalacia;
- Brittle bone disease;
- Communitated fractures (burst-fracture) and compression fractures;
- Tumours;
- Grade IV spondylolisthesis;
- Extensive epidural fibrosis;
- Active infectious process or major risk of infection;
- Weakened immune system;
- Fever or leukocytosis;
- Morbid obesity;
- Pregnancy;
- Hyperactivity;
- Mental illness;
- Allergy or intolerance to the implant materials;
- Any case not described in the indications.

7. Use

- Consult the surgical technique before use.

Medical device	Surgical approach
SHARK	Posterior approach and more specifically the transforaminal approach (TLIF)
PDP	Posterior approach and more specifically the unilateral approach (UNILIF)
PM CAGE®	Anterior approach

8. User profile

The implantation of this type of implant involves a standard surgical technique that has been used in spinal surgery for decades, but the surgeon must be duly trained in spinal surgery.

Reading the documentation and surgical techniques allows the surgeon to assimilate and understand the surgical risks associated with placement of the implants.

Sterilisation staff must refer to these instructions for use in particular to understand the cleaning and sterilisation procedures to be applied.

9. **Warning/Surgical precautions**

- The surgeon must be completely familiar with the implants, instruments and the surgical technique.
- The surgeon must take into account all of the factors likely to influence the performance of the device: level(s) to be instrumented, weight, activity level, concomitant disorders, etc.
- Although allergic reactions to titanium are very rare, it is recommended that patients with particularly strong allergies be checked beforehand for any allergies of this type.
- These implants must be used under perfectly sterile surgical conditions.
- Incorrect use of the instruments during implantation may lead to deterioration of the implants and therefore impact the expected performance.
- The instruments provided are designed to work perfectly with these devices. Using other instruments is therefore strongly contraindicated. This may incur risks for patients (see Sec. 10) and damage to the implants.
- Our implants are for single use only and must never be reused. Reuse may lead to reduced performance of the device, contamination and cross-infection.
- The implant must be disposed of if it has come into contact with blood, biological fluids or bodily tissues.
- The implant must be disposed of if it has been removed from the patient, damaged or if it has been used incorrectly.
- The implants are manufactured using Ti-6Al-4V titanium alloy and PEEK Optima™, which are non-magnetic materials and have a geometry that is unlikely to create significant induced currents. Moreover, as they are fixed to the bone and/or tissue, the implants are unlikely to be moved. Their safety, specifically in terms of heating and implant migration, has been assessed on the basis of bibliographical data comparing the available data on devices of similar composition, shape and use. This assessment concluded that they were safe for use in patients undergoing MRI examinations between 1.5 and 3 Tesla.
- As a precaution, it is advisable to inform the MRI operator of the presence of this implant before these examinations.
- Titanium implants can generate metal-related artefacts, which must be taken into account when performing and interpreting imaging examinations.
- The surgeon must be notified of any change (appearance, pain, etc.) at the site of the implant.
- Any type of accident must be reported to the surgeon, even if there are no visible exterior indications at the site of the implant.

10. **Adverse effects/complications**

The potential adverse effects or complications are:

- Severe pseudarthrosis;
- Superficial or deep infection at the site of the implant;
- Haematoma;
- Disorder of an adjacent disc;
- Neurological and/or vascular damage related to the surgical procedure;
- Displacement or expulsion of the implant requiring a repeat procedure;
- Fatigue fracture of the material;
- Death.

11. **Storage conditions for implants**

The implants must be stored in a dry place at room temperature away from sunlight and UV rays in their original packaging.

12. **Sterilisation**

The implants are delivered sterile under double vacuum packaging in a cardboard box, filmed and labelled. Sterilisation is performed by cobalt-60 irradiation of at least 25 kGy.

The sterilisation is valid for 5 years.

Caution: Do not use these implants in the following case: absence of vacuum, which proves a breach of the packaging, meaning loss of sterility.

The integrity of all levels of packaging must be checked prior to use along with the expiry date.

Additional labels are placed inside the packaging to ensure traceability in health institutions.

13. **Disposing of implants**

The removal and disposal of implants follow the same procedures used in the health institution as for hospital waste.

14. **Medical device vigilance**

Any complaint or notification of malfunction must be made immediately either by telephone, fax, email or post to NEURO FRANCE Implants.

NEURO FRANCE Implants and the competent authority of the Member State where the user and/or the patient is based must be notified of any serious incident involving the device.

15. **Lifetime**

The lifetime of these medical devices is 15 years. This consists of:

- a shelf life of 5 years;
- and an implantation period of 10 years.

Please note that this lifetime is based on mechanical, biological and clinical validations carried out on this medical device.

In reality, it is a permanent implantable device and is therefore not intended to be removed from the patient once implanted.

16. **Patient information**

The surgeon must warn the patient of the risks of surgery and inform them of the potential adverse effects.

The patient must be informed that the devices cannot reproduce the flexibility, strength, reliability or durability of normal healthy bone. Therefore, if the patient's profession or activities subject the implants to excessive stresses (e.g. lifting loads or making major muscular efforts, exposure to excessive mechanical vibration movements, etc.), the implants may break or be damaged, which may require another operation. The surgeon must therefore discuss all of the physical and psychological limitations inherent in the use of these devices with the patient. The surgeon must also explain the need for regular medical follow-up during the postoperative phase.

Smoking has been shown to impede patients' recovery after surgery and affect bone healing and repair. Surgeons must inform patients of this fact and warn them of the potential consequences.

Patients should be aware that they must declare that they have an implant before any exposure to electromagnetic and magnetic environments.

For further information, a patient information letter is available on the NEURO FRANCE Implants website.

Health institutions and/or practitioners are asked to complete the implant card provided by NEURO FRANCE Implants. This card is intended to inform the patient of the type and batch number of the implanted medical device, the health institution and the surgeon.

17. **Additional information**

Medical device	Date placed on the market
SHARK PEEK	July 2015
PDP	February 2015
PM CAGE® PEEK	May 2015

Information on the clinical results of these medical devices can be found in the “Summary of Safety and Clinical Performance”.
This report is available in the European Medical Device Database (Eudamed)*, where it is linked to the basic IUD-ID of the various interbody cages:

Medical device	Basic UDI-DI
SHARK PEEK	37005271CAGSHA03B2
PDP	37005271CAGPDP03BG
PM CAGE® PEEK	37005271CAGPMC03BE

The public Eudamed website is accessible via the following link: <https://ec.europa.eu/tools/eudamed>

Further information can be obtained from NEURO FRANCE Implants.

Explanation of the symbols used on the product labelling:

	Catalogue number		Manufacturer
	Batch code		Date of manufacture
	Do not re-use		Consult instructions for use
	Sterilised by irradiation		CE marking with the identification number of the notified body
	Use by date		Medical device
	Do not use if the packaging is damaged		Double sterile barrier system with protective external packaging

These implants are Class III medical devices and bear the CE marking. Our notified body is DQS

** Once the notified body has validated the report and the Eudamed site is fully operational.*

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