FIXATION RACHIDIENNE *SPINAL FIXATION*

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SPINAL FIXATION GENERIC GROUP INSTRUCTIONS FOR USE

Manufacturer: NEURO FRANCE Implants (NFI)
ZA LE BOURG – 25 RUE DES ÉCOLES
41160 LA VILLE AUX CLERCS – FRANCE
Telephone: (33) 02 54 80 90 90
Fax: (33) 02 54 80 83 33
Email: sav@neuro-france.net

Important information for the practitioner:

The practitioner must carefully read the instructions and guidelines, as well as the recommendations (presented in the surgical technique and in these instructions) before using the spinal fixation devices.

1. Description of the materials

- These instructions for use cover the spinal fixation generic group. It includes the dorso-lumbar fixation systems (KM, CTJ, G2S®, SES PLUS) and the cervical fixation devices (KM-C, JOC).

Medical device		Intended purpose
JOC	Occipito rod-plate Cervical screw	Cervico-cranial junction Cranial anchoring
KM-C	Polyaxial and polyaxial revision screw KM-C hook Pedicular hook Rod Cross link	Pedicular anchoring Laminar anchoring Pedicular anchoring Junction between two or more anchoring implants Junction between two rods
CTJ	Monoaxial screw Polyaxial screw Hook Pedicular hook Rod Link and connection system	Pedicular anchoring Pedicular anchoring Laminar anchoring Pedicular anchoring Junction between two or more anchoring implants Junction between two rods
KM	Monoaxial and monoaxial spondylo screw Polyaxial and polyaxial spondylo screw Hook Rod Link and connection system	Pedicular anchoring Pedicular anchoring Laminar anchoring Junction between two or more anchoring implants Junction between two rods
G2S®	Monoaxial thoracic, lumbar revision screw Polyaxial and polyaxial spondylo screw Cannulated screw Polyaxial thoracic, lumbar and spondylo cannulated screw Polyaxial thoraco-lumbar cannulated screw Hook Pedicular hook Precurved rod Ilio-sacral connector Link and connection system	Pedicular anchoring Pedicular anchoring Pedicular anchoring Pedicular anchoring Pedicular anchoring Laminar anchoring Pedicular anchoring Junction between two or more anchoring implants Junction between the rod and the pedicular screw Junction between two rods
SES PLUS	Monoaxial spondylo screw Ilio-sacral screw Sacral plate 3D washer and SES PLUS nut SES PLUS hook Straight or precurved rod Ø 5.5 mm Articulated cross link Connector	Pedicular anchoring Ilio-sacral anchoring Sacral screw-rod junction Locking the rod in the screw head Placing the rod onto the screw Junction between two or more anchoring implants Junction between two rods Junction between two rods

- The G2S® system can be used in combination with SpineGuard® DSG® technology.
- The JOC device is used in combination with the KM-C or CTJ systems.
- The table below presents the different assemblies possible depending on the devices used.

	Occipito rod-plate Ø3	Occipito rod-plate Ø5	Rod Ø3	Rod Ø5	Rod Ø5.5	Rod Ø6	Ilio-sacral screw	Nut
STABILAR screw	x	x						
KM-C screw			x					550125
CTJ screw				x				550300
KM screw				x	x	x		550300
G2S screw				x	x			550300
SES PLUS screw with SES PLUS hook					x			551004 (+530000)
SES PLUS sacral plate							x	551004 (+530000)
Cross link Ø3			x					550125
Articulated cross link				x	x	x		550250 550280
Cylindrical connector			x	x	x	x		550010 550011
Square connector			x	x	x	x		550010 550011
Ilio-sacral connector				x	x	x		550010 550011
KM-C hook			x					550125
CTJ hook				x				550300
KM hook				x	x	x		550300
G2S hook				x	x			550300

- These devices are made from Ti-6Al-4V titanium alloy, complying with ISO standard 5832-3.
- These devices are available in non-sterile condition.
- The batch number is marked on the implant. It is also printed on the product packaging.
- These implants are not intended to transport or introduce any substance into the human body.
- These implants do not contain any medicinal product, derivative of human blood or plasma, or tissue or cells of human or animal origin.
- Appropriate instruments are provided for implanting these devices. These ancillary products are reusable surgical instruments.
- These instruments are made of stainless steel or titanium. The handles are made of stainless steel, except for some which are made of silicone.
- They are delivered non-sterile. *Specific instructions for use are available for these instruments.*

2. Intended purpose

These osteosynthesis devices are intended to provide fixation between two or more vertebrae.

Medical device	Intended purpose
JOC	Posterior occipito-cervico-thoracic spinal arthrodesis
KM-C	Posterior cervical spinal arthrodesis
CTJ	Posterior cervico-thoracic spinal arthrodesis
KM	Posterior dorso-lumbar-sacral spinal arthrodesis
G2S	Posterior cervico-dorso-lumbar-sacral spinal arthrodesis
SES PLUS	Posterior or anterolateral dorso-lumbar-sacral spinal arthrodesis

The implants are not designed to withstand anatomical mechanical stresses for more than one year without bone fusion.

3. Expected clinical performance

- Stabilisation of the operated level(s) of the spine in over 90% of cases.
- In the case of scoliosis: improvement of the Cobb angle by at least 40%.

4. Expected clinical benefits

- Reduction in current clinical pain score (VAS) of at least 60% at 12 months postoperatively.
- Reduction in current clinical disability score (ODI) of at least 40% at 12 months postoperatively.

5. Indications

Medical device	Indications	Patient profile	Levels
JOC	Degenerative disc disease	Any adult patient (male or female) with one of the indications previously described and with no contraindications	Occiput - T3
KM-C	Fracture		C1 - C7
CTJ	Degenerative disc disease Scoliosis Fracture		C6 - T12
KM	Degenerative disc disease		T1 - S2
G2S	Spondylolisthesis		C6 - S2
SES PLUS	Scoliosis Fracture		T1 - S2

6. Contraindications

The following conditions may compromise the implantation (the list is not exhaustive):

- Severe osteoporosis
- Osteomalacia
- Brittle bone disease
- Comminuted fractures (burst fracture) and compression fractures
- Tumours
- Extensive epidural fibrosis
- Active infectious process or major risk of infection
- Weakened immune system
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Hyperactivity
- Mental illness
- Allergy or intolerance to the implant materials
- Any case not described in the indications

7. Use

- Consult the surgical technique before use.
- The choice of screw size must take into account the patient's body size and the level at which osteosynthesis is to be performed, in order to use pedicular implants of appropriate diameter and length.

Medical device	Surgical approach
JOC	Posterior approach
KM-C	
CTJ	
KM	
G2S	Classic posterior approach, percutaneous posterior approach and MIS (Minimally Invasive Surgery)
SES PLUS	Posterior or anterolateral approach

- When the G2S® system is used in combination with the DSG® technology by SpineGuard, refer to the SpineGuard surgical technique.
- When using fenestrated screws, Teknimed F20® cement or equivalent is recommended. It is a medium viscosity radiopaque surgical cement with a long working time. Please refer to the "Operator Information Sheet" for F20® or any other cement for further details.

8. User profile

The implantation of this type of implant involves a standard surgical technique that has been used in spinal surgery for decades, but the surgeon must be duly trained in spinal surgery. Additional training is necessary for the anterolateral approach.

Reading the documentation and surgical techniques allows the surgeon to assimilate and understand the surgical risks associated with placement of the implants.

Sterilisation staff must refer to these instructions for use in particular to understand the cleaning and sterilisation procedures to be applied.

9. Warning/Surgical precautions

- The surgeon must be completely familiar with the implants, instruments and the surgical technique.
- The surgeon must take into account all of the factors likely to influence the performance of the device: level(s) to be instrumented, weight, activity level, concomitant disorders, etc.
- Although allergic reactions to titanium are very rare, it is recommended that patients with particularly strong allergies be checked beforehand for any allergies of this type.
- The KM, CTJ, G2S®, KM-C, JOC and SES PLUS implants must be used under perfectly sterile surgical conditions.
- The components of the KM, CTJ, G2S®, KM-C, JOC and SES PLUS implants must not be used under any circumstances in combination with components other than those mentioned in paragraph 1. *Description of the materials*.
- Incorrect use of the instruments during implantation may lead to deterioration of the implants and therefore impact the expected performance.
- The instruments provided are designed to work perfectly with these devices. Using other instruments is therefore strongly contraindicated. This may incur risks for patients (see Sec. 10) and damage to the implants.
- Our implants are for single use only and must never be reused. Reuse may lead to reduced performance of the device, contamination and cross-infection.
- The implant must be disposed of if it has come into contact with blood, biological fluids or bodily tissues.
- The implant must be disposed of if it has been removed from the patient, damaged or if it has been used incorrectly.
- The implants are manufactured using Ti-6Al-4V titanium alloy, which is a non-magnetic material, and have a geometry that is unlikely to create significant induced currents. Moreover, as they are fixed to the bone and/or tissue, the implants are unlikely to be moved. Their safety, specifically in terms of heating and implant migration, has been assessed on the basis of bibliographical data comparing the available data on devices of similar composition, shape and use. This assessment concluded that they were safe for use in patients undergoing MRI examinations between 1.5 and 3 Tesla.
- As a precaution, it is advisable to inform the MRI operator of the presence of this implant before these examinations.
- Titanium implants can generate metal-related artefacts, which must be taken into account when performing and interpreting imaging examinations.
- The surgeon must be notified of any change (appearance, pain, etc.) at the site of the implant.
- Any type of accident must be reported to the surgeon, even if there are no visible exterior indications at the site of the implant.

10. Adverse effects/complications

The potential adverse effects or complications are:

- Neurological pain
- Durotomy
- Superficial or deep infection at the site of the implant
- Haematoma
- Disorder of an adjacent disc
- Neurological and/or vascular damage related to the surgical procedure
- Pedicular and/or laminar fracture
- Fatigue fracture of the material and/or migration of the material
- Death

11. Storage conditions for implants

The implants must be stored in a dry place at room temperature.

12. Cleaning

The implants must be cleaned and disinfected before being sterilised. The purpose of this step is to reduce the number of viable micro-organisms on the implants prior to sterilisation.

The chemicals used during cleaning and disinfection must be compatible with the implant. The implants must not come into contact with products containing chlorine, phosphorous, formalin, fluorine, bleaching agents, fatty-acid detergents or any strong acid or alkaline solutions. To guarantee effective cleaning and disinfection, please follow the cleaning product instructions.

The following method must be used for decontamination and cleaning:

Equipment required: Washer-disinfector, enzymatic cleaner/detergent with a pH between 8 and 10, such as Dentalzym5.

Programme validated with the SMEG WD4060D washer-disinfector.

1. Remove the implants from their packaging and place them in a basket on their respective racks.
2. Place the basket in the washer-disinfector using the validated process so as to avoid dead zones (Caution: do not stack multiple baskets on top of each other. Cleaning in a washer-disinfector must generally be carried out without a cover).
3. Follow the programme below:

SMEG WD4060D washer	Duration (in minutes)	Temperature	Cleaning agent	Quantity
Pre-wash	3	Unheated	Tap water	-
Wash cycle 1	5	35°C	Dentalzym5	6 mL/litre of water (Volume of water: 6 L)
Rinse cycle	-	Unheated	Tap water	-
Wash cycle 2	5	75°C	Dentalzym5	6 mL/litre of water (Volume of water: 6 L)
Rinse cycle	-	Unheated	Tap water	-
Neutralisation	3	-	Dental AC5	10 mL/litre of water (Volume of water: 9 L)
Rinse cycle	3	-	Tap water	-
Thermal disinfection	5	90°C	Reverse osmosis water	-
Drying cycle	25	110°C	-	-

4. When the programme ends, remove the basket from the washer-disinfector and visually confirm that the wash has been completed.
5. Place the basket in folded packaging (envelope or Pasteur).

13. Sterilisation

As the implants are delivered non-sterile, the health institution is responsible for sterilising them.

Specialised containers are available to staff for storing the implants.

Sterilisation must be carried out in an autoclave, using one of the following protocols.

Method 1	Method 2
Steam sterilisation (moist heat) Recommended time: 18 minutes (under saturated steam pressure) Recommended temperature: 134°C	Steam sterilisation (moist heat) Recommended time: 4 minutes (under saturated steam pressure) Recommended temperature: 132°C

It is recommended that operating theatre staff check the integrity of packaging after sterilisation in the autoclave.

14. Disposing of implants

The removal and disposal of implants follow the same procedures used in the health institution as for hospital waste.

15. Medical device vigilance

Any complaint or notification of malfunction must be made immediately either by telephone, fax, email or post to NEURO FRANCE Implants.

NEURO FRANCE Implants and the competent authority of the Member State where the user and/or the patient is based must be notified of any serious incident involving the device.

16. Lifetime

The lifetime of these medical devices is 20 years. This consists of:

- A maximum shelf life of 10 years
- And an implantation period of 10 years

Please note that this lifetime is based on mechanical, biological and clinical validations carried out on this medical device.

In reality, it is a permanent implantable device and is therefore not intended to be removed from the patient once implanted.

17. Patient information

The surgeon must warn the patient of the risks of surgery and inform them of the potential adverse effects.

The patient must be informed that the devices cannot reproduce the flexibility, strength, reliability or durability of normal healthy bone. Therefore, if the patient's profession or activities subject the implants to excessive stresses (e.g. lifting loads or making major muscular efforts, exposure to excessive mechanical vibration movements, etc.), the implants may break or be damaged, which may require another operation. The surgeon must therefore discuss all of the physical and psychological limitations inherent in the use of these devices with the patient. The surgeon must also explain the need for regular medical follow-up during the postoperative phase.

Smoking has been shown to impede patients' recovery after surgery and affect bone healing and repair. Surgeons must inform patients of this fact and warn them of the potential consequences.

Patients should be aware that they must declare that they have an implant before any exposure to electromagnetic and magnetic environments.

For further information, a patient information letter is available on the NEURO FRANCE Implants website.

Health institutions and/or practitioners are asked to complete the implant card provided by NEURO FRANCE Implants. This card is intended to inform the patient of the type and batch number of the implanted medical device, the health institution and the surgeon.

18. Additional information

Medical device	Date placed on the market
JOC	November 2006
KM-C	July 2006
KM	March 2005
CTJ	April 2010
G2S®	October 2014
SES PLUS	March 2019

Information on the clinical results of these medical devices can be found in the “Summary of Safety and Clinical Performance”. This report is available in the European Medical Device Database (Eudamed)*, where it is linked to the Basic UDI-DI of different spinal fixation devices:

Medical device	Basic UDI-DI
STABILAR® cervical screw	37005271VISCRA01TZ
JOC occipito rod-plate	37005271PBAJOC00FE
KM-C	37005271KMCSYSTEMG5
CTJ	37005271CTJSYSTEMJ9
KM	37005271KMLSYSTEMLW
G2S®	37005271G2SSYSTEM6B
SES PLUS	37005271SESSYSTEMRM

The public Eudamed website is accessible via the following link: <https://ec.europa.eu/tools/eudamed>

Further information can be obtained from NEURO FRANCE Implants.

Explanation of the symbols used on the product labelling:

	Catalogue number		Protective packaging
	Batch code		Manufacturer and date of manufacture
	Do not re-use		Consult instructions for use
	Non-sterile		CE marking with the identification number of the notified body
	Medical device		Unique Device Identifier

These implants are Class III medical devices and bear the CE marking. Our notified body is DQS.

**Once the notified body has validated the report and the Eudamed site is fully operational.*

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