CAGES NON-STERILE

En - LUMBAR & CERVICAL INTERBODY CAGES (PEEK and TITANE, NON-STERILE)
INSTRUCTIONS FOR USEp.2

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neurofrance@neuro-france.net

www.neuro-france.net

LUMBAR & CERVICAL INTERBODY CAGES (PEEK and TITANE, NON-STERILE) INSTRUCTIONS FOR USE

Manufacturer: NEURO FRANCE Implants (NFI)
ZA LE BOURG – 25 RUE DES ECOLES
41160 LA VILLE AUX CLERCS – FRANCE
Telephone: (33) 02 54 80 90 90
Fax: (33) 02 54 80 83 33
Email: sav@neuro-france.net

Important information for the practitioner:

The practitioner must carefully read the instructions and guidelines, as well as the recommendations (presented in the surgical technique and in these instructions) before using a lumbar or cervical interbody cage.

1. Description of the materials

- These instructions for use cover the interbody cages generic group. It includes lumbar interbody cages (SHARK, PDP, AVENIR, BANANA, EUROPA®) and cervical interbody cages (PM CAGE®).

Medical device	Material
SHARK	PEEK Optima™ (with radiological markers made of Ti-6Al-4V titanium alloy)
	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
BANANA	PEEK Optima™ (with radiological markers made of Ti-6Al-4V titanium alloy)
	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
PDP	PEEK Optima™ (with radiological markers made of Ti-6Al-4V titanium alloy)
AVENIR	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
EUROPA®	Ti-6Al-4V titanium alloy (standard ISO 5832-3)
PM CAGE®	PEEK Optima™ (with radiological markers made of Ti-6Al-4V titanium alloy)
	Ti-6Al-4V titanium alloy (standard ISO 5832-3)

- These cages are available in non-sterile condition.
- The batch number is marked on the implant. It is also printed on the product packaging.
- These implants are not intended to transport or introduce any substance into the human body.
- These implants do not contain any medicinal product, derivative of human blood or plasma, or tissue or cells of human or animal origin.
- Appropriate instruments are provided for implanting these devices. These ancillary products are reusable surgical instruments.
- These instruments are made of stainless steel or titanium. The handles are made of stainless steel, except for some which are made of silicone.
- They are delivered non-sterile. *Specific instructions for use are available for these instruments.*

2. Intended purpose

These cages are designed for interbody fusion. These implants restore the height of the intervertebral space.

Note that lumbar interbody cages were not designed as "independent" implants, so they require the use of posterior instrumentation (such as pedicle screws or plates).

3. Expected clinical performance

- Fusion rate greater than 85% at 12 months postoperatively.
- Increase in disc height of at least 15%.

4. Expected clinical benefits

- Reduction in current clinical pain score (VAS) of at least 60% at 12 months postoperatively.
- Reduction in current clinical disability score (NDI/ODI) of at least 40% at 12 months postoperatively.

5. Indications

Medical device	Indications	Patient profile	Levels
SHARK	Degenerative disc disease Spondylolisthesis	Any adult patient (male or female) with one of the indications previously described and no contraindications.	L1 - S1
BANANA			
PDP			
AVENIR			
EUROPA®			
PM CAGE®	Degenerative disc disease Cervico-brachial neuralgia (myelopathy - radiculopathy)		C2 – T1

6. Contraindications

The following conditions may compromise the implantation (the list is not exhaustive):

- Severe osteoporosis;
- Osteomalacia;
- Brittle bone disease;
- Communitated fractures (burst-fracture) and compression fractures;
- Tumours;
- Grade IV spondylolisthesis;
- Extensive epidural fibrosis;
- Active infectious process or major risk of infection;
- Weakened immune system;
- Fever or leukocytosis;
- Morbid obesity;
- Pregnancy;
- Hyperactivity;
- Mental illness;
- Allergy or intolerance to the implant materials;
- Any case not described in the indications.

7. Use

- Consult the surgical technique before use.

Medical device	Surgical approach
SHARK	Posterior approach and more specifically the transforaminal approach (TLIF)
EUROPA®	
BANANA	Posterior approach and more specifically the transforaminal or unilateral approach
PDP	Posterior approach and more specifically the unilateral approach (UNILIF)
AVENIR	
PM CAGE®	Anterior approach

8. User profile

The implantation of this type of implant involves a standard surgical technique that has been used in spinal surgery for decades, but the surgeon must be duly trained in spinal surgery.

Reading the documentation and surgical techniques allows the surgeon to assimilate and understand the surgical risks associated with placement of the implants.

Sterilisation staff must refer to these instructions for use in particular to understand the cleaning and sterilisation procedures to be applied.

9. Warning/Surgical precautions

- The surgeon must be completely familiar with the implants, instruments and the surgical technique.
- The surgeon must take into account all of the factors likely to influence the performance of the device: level(s) to be instrumented, weight, activity level, concomitant disorders, etc.
- Although allergic reactions to titanium are very rare, it is recommended that patients with particularly strong allergies be checked beforehand for any allergies of this type.
- These implants must be used under perfectly sterile surgical conditions.
- Incorrect use of the instruments during implantation may lead to deterioration of the implants and therefore impact the expected performance.
- The instruments provided are designed to work perfectly with these devices. Using other instruments is therefore strongly contraindicated. This may incur risks for patients (see Sec. 10) and damage to the implants.
- Our implants are for single use only and must never be reused. Reuse may lead to reduced performance of the device, contamination and cross-infection.
- The implant must be disposed of if it has come into contact with blood, biological fluids or bodily tissues.
- The implant must be disposed of if it has been removed from the patient, damaged or if it has been used incorrectly.
- The implants are manufactured using Ti-6Al-4V titanium alloy and PEEK Optima™, which are non-magnetic materials and have a geometry that is unlikely to create significant induced currents. Moreover, as they are fixed to the bone and/or tissue, the implants are unlikely to be moved. Their safety, specifically in terms of heating and implant migration, has been assessed on the basis of bibliographical data comparing the available data on devices of similar composition, shape and use. This assessment concluded that they were safe for use in patients undergoing MRI examinations between 1.5 and 3 Tesla.
- As a precaution, it is advisable to inform the MRI operator of the presence of this implant before these examinations.
- Titanium implants can generate metal-related artefacts, which must be taken into account when performing and interpreting imaging examinations.
- The surgeon must be notified of any change (appearance, pain, etc.) at the site of the implant.
- Any type of accident must be reported to the surgeon, even if there are no visible exterior indications at the site of the implant.

10. Adverse effects/complications

The potential adverse effects or complications are:

- Severe pseudarthrosis;
- Superficial or deep infection at the site of the implant;
- Haematoma;
- Disorder of an adjacent disc;
- Neurological and/or vascular damage related to the surgical procedure;
- Displacement or expulsion of the implant requiring a repeat procedure;
- Fatigue fracture of the material;
- Death.

11. Storage conditions for implants

The implants must be stored in a dry place at room temperature.

12. Cleaning

The implants must be cleaned and disinfected before being sterilised. The purpose of this step is to reduce the number of viable micro-organisms on the implants prior to sterilisation.

The chemicals used during cleaning and disinfection must be compatible with the implant. The implants must not come into contact with products containing chlorine, phosphorous, formalin, fluorine, bleaching agents, fatty-acid detergents or any strong acid or alkaline solutions.

To guarantee effective cleaning and disinfection, please follow the cleaning product instructions.

The following method must be used for decontamination and cleaning:

Equipment required: Washer-disinfector, enzymatic cleaner/detergent with a pH between 8 and 10, such as Dentalzym5.

Programme validated with the SMEG WD4060D washer-disinfector.

- Remove the implants from their packaging and place them in a basket on their respective racks.
- Place the basket in the washer-disinfector using the validated process so as to avoid dead zones (Caution: do not stack multiple baskets on top of each other. Cleaning in a washer-disinfector must generally be carried out without a cover).
- Follow the programme below:

SMEG WD4060D washer	Duration (in minutes)	Temperature	Cleaning agent	Quantity
Pre-wash	3	Unheated	Tap water	-
Wash cycle 1	5	35°C	Dentalzym5	6mL/litre of water (Volume of water: 6L)
Rinse cycle	-	Unheated	Tap water	-
Wash cycle 2	5	75°C	Dentalzym5	6mL/litre of water (Volume of water: 6L)
Rinse cycle	-	Unheated	Tap water	-
Neutralisation	3	-	Dental AC5	10mL/litre of water (Volume of water: 9L)
Rinse cycle	3	-	Tap water	-
Thermal disinfection	5	90°C	Reverse osmosis water	-
Drying cycle	25	110°C	-	-

4. When the programme ends, remove the basket from the washer-disinfector and visually confirm that the wash has been completed.
5. Place the basket in folded packaging (envelope or Pasteur).

13. Sterilisation

As the implants are delivered non-sterile, the health institution is responsible for sterilising them.

Specialised containers are available to staff for storing the implants.

Sterilisation must be carried out in an autoclave, using one of the following protocols

Method 1	Method 2
Steam sterilisation (moist heat) Recommended time: 18 minutes (under saturated steam pressure) Recommended temperature: 134°C	Steam sterilisation (moist heat) Recommended time: 4 minutes (under saturated steam pressure) Recommended temperature: 132°C

It is recommended that operating theatre staff check the integrity of the packaging after sterilisation in the autoclave.

14. Disposing of implants

The removal and disposal of implants follow the same procedures used in the health institution as for hospital waste.

15. Medical device vigilance

Any complaint or notification of malfunction must be made immediately either by telephone, fax, email or post to NEURO FRANCE Implants.

NEURO FRANCE Implants and the competent authority of the Member State where the user and/or the patient is based must be notified of any serious incident involving the device.

16. Lifetime

The lifetime of these medical devices is 20 years. This consists of:

- a maximum shelf life of 10 years;
- and an implantation period of 10 years.

Please note that this lifetime is based on mechanical, biological and clinical validations carried out on this medical device.

In reality, it is a permanent implantable device and is therefore not intended to be removed from the patient once implanted.

17. Patient information

The surgeon must warn the patient of the risks of surgery and inform them of the potential adverse effects.

The patient must be informed that the devices cannot reproduce the flexibility, strength, reliability or durability of normal healthy bone. Therefore, if the patient's profession or activities subject the implants to excessive stresses (e.g. lifting loads or making major muscular efforts, exposure to excessive mechanical vibration movements, etc.), the implants may break or be damaged, which may require another operation. The surgeon must therefore discuss all of the physical and psychological limitations inherent in the use of these devices with the patient. The surgeon must also explain the need for regular medical follow-up during the postoperative phase.

Smoking has been shown to impede patients' recovery after surgery and affect bone healing and repair. Surgeons must inform patients of this fact and warn them of the potential consequences.

Patients should be aware that they must declare that they have an implant before any exposure to electromagnetic and magnetic environments.

For further information, a patient information letter is available on the NEURO FRANCE Implants website.

Health institutions and/or practitioners are asked to complete the implant card provided by NEURO FRANCE Implants. This card is intended to inform the patient of the type and batch number of the implanted medical device, the health institution and the surgeon.

18. Additional information

Medical device	Date placed on the market
SHARK PEEK	April 2010
SHARK TITANE	February 2019
BANANA PEEK	June 2017
BANANA TITANE	February 2019
PDP	June 2012
AVENIR	February 2019
EUROPA®	September 2007
PM CAGE® PEEK	March 2005
PM CAGE® TITANE	February 2017

Information on the clinical results of these medical devices can be found in the "Summary of Safety and Clinical Performance". This report is available in the European Medical Device Database (Eudamed)*, where it is linked to the basic IUD-ID of the various interbody cages:

Medical device	Basic UDI-DI
SHARK PEEK	37005271CAGSHA01AW
SHARK TITANE	37005271CAGSHA02AY
BANANA PEEK	37005271CAGBAN015K
BANANA TITANE	37005271CAGBAN025M
PDP	37005271CAGPDP01BC
AVENIR	37005271CAGAVE028G
EUROPA®	37005271CAGEUR02BN
PM CAGE® PEEK	37005271CAGPMC01BA
PM CAGE® TITANE	37005271CAGPMC02BC

The public Eudamed website is accessible via the following link: <https://ec.europa.eu/tools/eudamed>

Further information can be obtained from NEURO FRANCE Implants.

Explanation of the symbols used on the product labelling:

	Catalogue number		Protective packaging
	Batch code		Manufacturer and date of manufacture
	Do not re-use		Consult instructions for use
	Non-sterile		CE marking with the identification number of the notified body
	Medical device		Unique Device Identifier

These implants are Class III medical devices and bear the CE marking. Our notified body is DQS

* Once the notified body has validated the report and the Eudamed site is fully operational.

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