G2S® STERILE

De - POSTERIORES DORSAL-LUMBALES OSTEOSYNTHESYSTEM G2S® Global Spine System
GEBRAUCHSANWEISUNGp.2

En - POSTERIOR THORACO-LUMBAR OSTEOSYNTHESIS SYSTEM G2S® Global Spine System
INSTRUCTIONS FOR USEp.5

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Important information for the practitioner:

The practitioner must carefully read the instructions and guidelines, as well as the recommendations (presented in the surgical technique and in these instructions) before using the G2S® system.

1. Description of the materials

- These instructions for use cover some of the components of the G2S® system which are available in sterile condition.
- These are the following components:

Medical device	Reference	Intended purpose
Cannulated polyaxial thoraco-lumbar screw with nut	2G305530S, 2G305535S, 2G305540S, 2G305545S, 2G305550S, 2G306535S, 2G306540S, 2G306545S, 2G306550S, 2G307530S, 2G307535S, 2G307540S, 2G307545S, 2G307550S	Pedicular anchoring
Cannulated polyaxial lumbar screw with nut	2G325530S, 2G325535S, 2G325540S, 2G325545S, 2G325550S, 2G326530S, 2G326535S, 2G326540S, 2G326545S, 2G326550S, 2G326555S, 2G327530S, 2G327535S, 2G327540S, 2G327545S, 2G327550S	
Precurved rod	4G525040S, 4G525045S, 4G525055S, 4G525065S, 4G525075S, 4G525090S, 4G525105S, 4G525125S	Junction between two or more anchoring implants
Precurved percutaneous rod	411505S, 411506S, 411507S, 411508S, 411511S, 411515S, 411516S, 411517S	

- These devices are made from Ti-6Al-4V titanium alloy, complying with ISO standard 5832-3.
- The batch number is marked on the implant. It is also printed on the product packaging.
- These implants are not intended to transport or introduce any substance into the human body.
- These implants do not contain any medicinal product, derivative of human blood or plasma, or tissue or cells of human or animal origin.
- Appropriate instruments are provided for implanting these devices. These ancillary products are reusable surgical instruments.
- These instruments are made of stainless steel or titanium. The handles are made of stainless steel, except for some which are made of silicone.
- They are delivered non-sterile. *Specific instructions for use are available for these instruments.*

2. Intended purpose

The G2S® system is intended to provide posterior thoraco-lumbar spinal arthrodesis by ensuring fixation between two or more vertebrae.

The implants are not designed to withstand anatomical mechanical stresses for more than one year without bone fusion.

3. Expected clinical performance

- Stabilisation of the operated level(s) of the spine in over 90% of cases.
- In the case of scoliosis: improvement of the Cobb angle by at least 40%.

4. Expected clinical benefits

- Reduction in current clinical pain score (VAS) of at least 60% at 12 months postoperatively.
- Reduction in current clinical disability score (ODI) of at least 40% at 12 months postoperatively.

5. Indications

Medical device	Indications	Patient profile	Levels
G2S®	Degenerative disc disease Spondylolisthesis Scoliosis Fracture	Any adult patient (male or female) with one of the previously described conditions and with no contraindications.	T10 – S1

6. Contraindications

The following conditions may compromise the implantation (the list is not exhaustive):

- Severe osteoporosis
- Osteomalacia
- Brittle bone disease
- Comminuted fractures (burst fracture) and compression fractures
- Tumours
- Extensive epidural fibrosis
- Active infectious process or major risk of infection
- Weakened immune system
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Hyperactivity
- Mental illness
- Allergy or intolerance to the implant materials
- Any case not described in the indications

7. Use

- Consult the surgical technique before use.
- The choice of screw size must take into account the patient's body size and the level at which osteosynthesis is to be performed, in order to use pedicular implants of appropriate diameter and length.
- The surgical approach is classic posterior approach, percutaneous posterior approach and MIS (Minimally Invasive Surgery).

8. User profile

The implantation of this type of implant involves a standard surgical technique that has been used in spinal surgery for decades, but the surgeon must be duly trained in spinal surgery.

Reading the documentation and surgical techniques allows the surgeon to assimilate and understand the surgical risks associated with placement of the implants.

Sterilisation staff must refer to these instructions for use in particular to understand the cleaning and sterilisation procedures to be applied.

9. Warning/Surgical precautions

- The surgeon must be completely familiar with the implants, instruments and the surgical technique.
- The surgeon must take into account all of the factors likely to influence the performance of the device: level(s) to be instrumented, weight, activity level, concomitant disorders, etc.
- Although allergic reactions to titanium are very rare, it is recommended that patients with particularly strong allergies be checked beforehand for any allergies of this type.
- The G2S® implants must be used under perfectly sterile surgical conditions.
- Incorrect use of the instruments during implantation may lead to deterioration of the implants and therefore impact the expected performance.
- The instruments provided are designed to work perfectly with these devices. Using other instruments is therefore strongly contraindicated. This may incur risks for patients (see Sec. 10) and damage to the implants.
- Our implants are for single use only and must never be reused. Reuse may lead to reduced performance of the device, contamination and cross-infection.
- The implant must be disposed of if it has come into contact with blood, biological fluids or bodily tissues.
- The implant must be disposed of if it has been removed from the patient, damaged or if it has been used incorrectly.
- The implants are manufactured using Ti-6Al-4V titanium alloy, which is a non-magnetic material, and have a geometry that is unlikely to create significant induced currents. Moreover, as they are fixed to the bone and/or tissue, the implants are unlikely to be moved. Their safety, specifically in terms of heating and implant migration, has been assessed on the basis of bibliographical data comparing the available data on devices of similar composition, shape and use. This assessment concluded that they were safe for use in patients undergoing MRI examinations between 1.5 and 3 Tesla.
- As a precaution, it is advisable to inform the MRI operator of the presence of this implant before these examinations.
- Titanium implants can generate metal-related artefacts, which must be taken into account when performing and interpreting imaging examinations.
- The surgeon must be notified of any change (appearance, pain, etc.) at the site of the implant.
- Any type of accident must be reported to the surgeon, even if there are no visible exterior indications at the site of the implant.

10. Adverse effects/complications

The potential adverse effects or complications are:

- Neurological pain
- Durotomy
- Superficial or deep infection at the site of the implant
- Haematoma
- Disorder of an adjacent disc
- Neurological and/or vascular damage related to the surgical procedure
- Pedicular and/or laminar fracture
- Fatigue fracture of the material and/or migration of the material
- Death

11. Storage conditions for implants

The implants must be stored in a dry place at room temperature away from sunlight and UV rays in their original packaging.

12. Sterilisation

The implants are delivered sterile under double vacuum packaging in a cardboard box, filmed and labelled. Sterilisation is performed by cobalt-60 irradiation of at least 25 kGy.

The sterilisation is valid for 5 years.

Caution: Do not use these implants in the following case: absence of vacuum, which proves a breach of the packaging, meaning loss of sterility.

The integrity of all levels of packaging must be checked prior to use along with the expiry date.

Additional labels are placed inside the packaging to ensure traceability in health institutions.

13. Disposing of implants

The removal and disposal of implants follow the same procedures used in the health institution as for hospital waste.

14. Medical device vigilance

Any complaint or notification of malfunction must be made immediately either by telephone, fax, email or post to NEURO FRANCE Implants.

NEURO FRANCE Implants and the competent authority of the Member State where the user and/or the patient is based must be notified of any serious incident involving the device.

15. Lifetime

The lifetime of these medical devices is 20 years. This consists of:

- A maximum shelf life of 10 years
- And an implantation period of 10 years

Please note that this lifetime is based on mechanical, biological and clinical validations carried out on this medical device. In reality, it is a permanent implantable device and is therefore not intended to be removed from the patient once implanted.

16. Patient information

The surgeon must warn the patient of the risks of surgery and inform them of the potential adverse effects.

The patient must be informed that the devices cannot reproduce the flexibility, strength, reliability or durability of normal healthy bone. Therefore, if the patient's profession or activities subject the implants to excessive stresses (e.g. lifting loads or making major muscular efforts, exposure to excessive mechanical vibration movements, etc.), the implants may break or be damaged, which may require another operation. The surgeon must therefore discuss all of the physical and psychological limitations inherent in the use of these devices with the patient. The surgeon must also explain the need for regular medical follow-up during the postoperative phase.

Smoking has been shown to impede patients' recovery after surgery and affect bone healing and repair. Surgeons must inform patients of this fact and warn them of the potential consequences.

Patients should be aware that they must declare that they have an implant before any exposure to electromagnetic and magnetic environments.

For further information, a patient information letter is available on the NEURO FRANCE Implants website.

Health institutions and/or practitioners are asked to complete the implant card provided by NEURO FRANCE Implants. This card is intended to inform the patient of the type and batch number of the implanted medical device, the health institution and the surgeon.

17. Additional information

Date placed on the market: October 2014

Information on the clinical results of these medical devices can be found in the "Summary of Safety and Clinical Performance". This report is available in the European Medical Device Database (Eudamed)*, where it is linked to the Basic UDI-DI of the system: 37005271G2SSYSTEM6B.

The public Eudamed website is accessible via the following link: <https://ec.europa.eu/tools/eudamed>

Further information can be obtained from NEURO FRANCE Implants.

Explanation of the symbols used on the product labelling:

	Catalogue number		Manufacturer
	Batch code		Date of manufacture
	Do not re-use		Consult instructions for use
	Sterilised by irradiation		CE marking with the identification number of the notified body
	Use by date		Medical device
	Do not use if the packaging is damaged		Double sterile barrier system with protective external packaging

These implants are Class III medical devices and bear the CE marking. Our notified body is DQS.

**Once the notified body has validated the report and the Eudamed site is fully operational.*

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Dispositif médical	Référence	Usage revendiqué
Vis canulée polyaxiale thoraco-lombaire avec écrou	2G305530S, 2G305535S, 2G305540S, 2G305545S, 2G305550S, 2G306535S, 2G306540S, 2G306545S, 2G306550S, 2G307530S, 2G307535S, 2G307540S, 2G307545S, 2G307550S	Ancrage pédiculaire
Vis canulée polyaxiale lombaire avec écrou	2G325530S, 2G325535S, 2G325540S, 2G325545S, 2G325550S, 2G326530S, 2G326535S, 2G326540S, 2G326545S, 2G326550S, 2G326555S, 2G327530S, 2G327535S, 2G327540S, 2G327545S, 2G327550S	
Barre pré-cambree	4G525040S, 4G525045S, 4G525055S, 4G525065S, 4G525075S, 4G525090S, 4G525105S, 4G525125S	Jonction entre deux ou plusieurs implants
Barre percutanée pré-cambree	411505S, 411506S, 411507S, 411508S, 411511S, 411515S, 411516S, 411517S	d'ancrage

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