



2021



+33 (0)2 54 80 90 90



Z.A. «Le Bourg» - 25 rue des Ecoles  
41160 La Ville aux Clercs - FRANCE



[www.neuro-france.net](http://www.neuro-france.net)






**TRIONYX®**

► TOTAL STERNAL PROSTHESIS

TRIONYX® PROSTHESIS..... 04



**ThoRib®**

► THORACIC  
OSTEOSYNTHESIS

ThoRib®..... 05

ThoRib® PECTUS..... 06



### ► Indications

- Sternal and/or manubrial bone tumor.

### ► Contents of the system

- TRIONYX® sternal plates (range of 6 sizes)
- Reinforced sternal staples
- Right precurved sternal staples
- Left precurved sternal staples
- Anti-vibration screws

### ► Specifications

After a sternal resection, the TRIONYX® prosthesis recreates the thoracic rigidity necessary to ensure normal respiratory mechanics and thus eliminate the risk of paradoxal breathing.

The TRIONYX® implant allows a stable and solid reconstruction of the sternum and/or manubrium after total or partial sternectomy.

TRIONYX® protects vital and mediastinal organs such as the heart, large vessels and lungs.





## THORACIC OSTEOSYNTHESIS SYSTEM



### ► Indications

- Severe thoracic trauma
- Thoracic reconstruction
- Rib stabilization after invasive breast tumor resection
- Malignant lung tumor necessitating rib resection
- Gyboplasty.

### ► Contents of the system

- Costal staples
- Large costal staples
- Simple costal staples
- Large simple costal staples
- Double staples (inter-axis from 30 to 100 mm)
- ThoRib® rods
- Anti-vibration screws
- ThoRib® connector

### ► Specifications

The THORIB® system solves the problems of thoracic osteosynthesis.

The **ThoRib® emergency kit** allows rapid management of chest trauma.

**The ThoRib® connector** associated with spinal fixation systems such as CTJ and G2S® makes it possible to protect the dura mater from pulmonary compression in the case of Pancoast tumor.

The ThoRib® system is available in sterile or non-sterile packaging.





### ► Indications

- Chest wall distortion and malformation treatment
- Pectus excavatum, carinatum and arcuatum

### ► Contents of the system

- PECTUS staples (thickness : 2,50 mm)
- PECTUS rods (thickness : from 2,50 mm to 3,50 mm)
- Anti-vibration screws


### ► Specifications

The ThoRib® PECTUS is a system of sternal rods designed to treat congenital deformities of the chest wall.

Thicker than a standard ThoRib® implant, ThoRib® PECTUS rods and staples have been specially developed to treat severe thoracic deformities. They allow a solid maintenance of the chest wall and a suppression of cardiac and pulmonary compression in the case of a PECTUS EXCAVATUM.

The PECTUS system also makes it possible to reduce the physical complexes of the patients related to these deformations.

The ThoRib® PECTUS system is available in sterile or non-sterile packaging.





Your partner in thoracic surgery

NEURO FRANCE Implants



Your partner in thoracic surgery

NEURO FRANCE Implants



Z.A «Le Bourg» - 25 rue des Ecoles  
41160 La Ville aux Clercs - FRANCE

Phone : +33 (0)2 54 80 90 90 / Fax : +33 (0)2 54 80 83 33  
neurofrance@neuro-france.net

[www.neuro-france.net](http://www.neuro-france.net)



CE  
459